

# Manganese


## 1. Determination of Manganese in Water

You have to determine the amount of manganese II in water. For this purpose, prepare a standard calibration curve with potassium permanganate. The aqueous manganese II salt solution is oxidized with ammonium peroxy disulfate, nitric acid and silver nitrate to give the Manganese VII salt (see reaction 1). The given water sample is measured three times.


The potassium permanganate solution, facing the fact of its ever changing concentration, is titrated with the primary standard (see reaction 2) and the exact substance amount concentration  $c(\text{KMnO}_4)$  is calculated.

## 2. Reaction Schemes:

### Reaction 1



### Reaction 2



## 3. Reagents for the Manganese Analysis:

Reagent	Amount	Information	Safety advice
Ammoniumperoxodisulfate	15g	Oxidation agent	H: 272-302-315-319-335-334-317 P:280-305+351+338-302+352-304+341-342+311
Nitric acid 65%	9*2.000mL	Oxidation agent	H: 272-314 P: 220-280-305+351+338-310
Silver nitrate 0.1M	9*2.000mL	Catalyst	272-314-410 P: 273-280-301+330+331-305+351+338-309+310
Sample with Mn II salt	3*10.00mL		H: 373-411 P: 273-314
<i>di</i> -Sodium oxalate	3*100mg	Primary standard	H: 302-312 P: 262

## Analytical Tasks

Muttenz 2015

Potassium permanganate standard solution 0.02M	6.5mL calib. 3*15mL		H: 272-302-314-410 P: 220-273-280-305+351+338-310-501
Sulfuric acid about 20% in water	3*20mL	Acidifying agent	H: 314-290 P: 280-301+330+331-309-310-305+351+338

## 4. Equipment for Manganese Analysis

Analytical balance	1 4all	Wash bottle with water <i>nanopure</i>	1
Lab balance	1 4all	2.00mL volumetric pipet	2
Bin or flask with 2.5L water <i>nanopure</i>	1 4all	4.00mL volumetric pipet	2
Automatic burette with KMnO <sub>4</sub> 0.02M	1 4all	Beaker 400mL	1
Piston stroke pipette 1000μL	1 4all	100mL Erlenmeyer flask with stopper	3
Piston stroke pipette 2500μL	1 4all	Magnetic stirrer with heat plate	1
Beaker 100mL	4	Magnetic stirring bar	6
Volumetric flask 100mL with stopper	9	Magnetic Teflon stick	1
Marker pen	1	Thermometer (20-100°C must be visible)	1
Spatula	3	PS (Polystyrol) cuvettes	9
Funnel for transferring prepared solutions	2	Cuvette rack	1
10.00mL volumetric pipet	1	Measuring cylinder 25mL	1
250mL Erlenmeyer flask	1		

Pasteur pipettes as required.

## 5. Procedure:

### Dissolved silver salts

Some of your volumetric flasks will contain silver nitrate salts. After use, they must be collected in special containers and never, under any circumstances, poured into the sink. They will be treated with sulfide later. (See number 13 *General Notes*)

### Preparation of Ammonium Peroxo Disulfate Solution

Weigh out 15g Ammoniumperoxodisulfate in a 250mL Erlenmeyer flask, dilute it with not more than 50mL water *nanopure* and pour it quantitatively into a 100.0mL measuring flask. Then the flask is filled with water *nanopure* to the mark.

### Preparation of the Samples

You have to prepare and measure the sample three times.

10.00mL of your sample are transferred by means of a volumetric pipet into a 100mL Erlenmeyer flask and diluted with approximately 10mL water. With a volumetric pipet transfer 2.00mL nitric acid  $w(\text{HNO}_3) = 0.65\text{g/g}$ , 4.00mL of the previous prepared ammonium peroxy disulfate solution. Then add 2'000 $\mu\text{L}$  silver nitrate  $c(\text{AgNO}_3) = 0.1\text{mol/L}$  with a 2'500 $\mu\text{L}$  piston stroke pipet (Note 1). Heat the solution up to at least 75°C (Note 2). The Erlenmeyer flask is stoppered and under occasional stirring or swiveling the violet mixture is cooled down to ambient temperature during circa 90 minutes. You can have all three Erlenmeyer flasks on the same stirrer plate for cooling. The cool solution is then transferred quantitatively and without the magnetic stirring bar into a 100mL volumetric flask. All volumetric flasks are filled with water *nanopure* to the mark (Note 3) and well homogenized.

### Preparation of Calibration Curve

Add 400 $\mu\text{L}$ , 850 $\mu\text{L}$ , 1'300 $\mu\text{L}$ , 1'750 $\mu\text{L}$  and 2'200 $\mu\text{L}$  of potassium permanganate  $c(\text{KMnO}_4) \approx 0.02\text{mol/L}$  into 100.0mL volumetric flasks, which are reasonably marked. Into each flask, transfer with a volumetric pipet 2.00mL nitric acid  $w(\text{HNO}_3) = 0.65\text{g/g}$ , 4.00mL of the previous prepared ammonium peroxy disulfate solution; Then add 2'000 $\mu\text{L}$  silver nitrate  $c(\text{AgNO}_3) = 0.1\text{mol/L}$  with a 2'500 $\mu\text{L}$  piston stroke pipet. All volumetric flasks are filled with water *nanopure* to the mark and well homogenized.

### Blank solution

Transfer with a volumetric pipet 2.00mL nitric acid  $w(\text{HNO}_3) = 0.65\text{g/g}$ , 4.00mL of the previous prepared ammonium peroxy disulfate solution into a 100mL volumetric flask. Then add 2'000 $\mu\text{L}$  silver nitrate  $c(\text{AgNO}_3) = 0.1\text{mol/L}$  with a 2'500 $\mu\text{L}$  piston stroke pipet. The volumetric flask is filled with water *nanopure* to the mark and well homogenized.

### Photometry and Results

After homogenizing the absorbances are measured photometrically at 526nm against the blank solution. The photospectrometer is already prepared for measuring at this wavelength. There is only one photospectrometer in the lab. In view of this fact, you have to prepare the PS cuvettes and measure the samples without delay, if it's your turn.

Note all absorbances found, as well as the appropriate concentration in a reasonable manner in your protocol.

### Calculations

Calculate the mass concentration  $\beta(\text{Mn})$  in your sample with the aid of your calibration curve. Involve the result of the volumetric analysis.

A linear equation could be helpful. Alternatively you can calculate the mass concentrations  $\beta(\text{Mn})$  in your samples with one point calibrations using close absorbances of the calibration curve. Also, you can use a blank excel-sheet on a computer or calculate with the aid of a hand calculator and millimeter paper.

## 6. Volumetric Analysis of Potassium Permanganate

Potassium permanganate standard solution is known for the fact that it tends to change its concentration slightly but steadily. In order to measure the exact amount of potassium permanganate we measure the exact substance amount concentration  $c(\text{KMnO}_4)$  with the aid of the primary standard di-sodium oxalate.

Weigh three times approximately 100mg di-natrium oxalate directly into a 100mL beaker. Note the exact amount of di-natrium oxalate primary standard. Add 20mL water *nanopure* and 20mL sulfuric acid 20%. Stir by means of a magnetic stirrer until a clear solution is present. This solution is then warmed on a hot magnetic stirring plate up to at least 50 to 80°C (Note 2).

The warm solution is titrated with the standard solution to the endpoint visible on the change from colorless to a slight pink color (Note 4).

Calculate the exact substance amount concentration of your standard solution.

## 7. Notes

**Note 1:** There are four 50 or 100mL beakers. You can fill them with the reagent solutions and use them for pipetting conveniently. Take potassium permanganate solution from the automatic burette. Never pour reagent solutions back into the storage bottle. Pour waste reagents in the silver waste or, if there is no silver present, under dilution with fresh water into the sink.

**Note 2:** There is no harm, if the mixture is heated up to a slight boiling. Avoid heavy boiling; you could bias your analysis if there is some splashing. If you use a thermometer, rinse it well with water *nanopure* into the solution.

**Note 3:** There is a 400mL beaker present. You can fill it with a suitable amount of water *nanopure*. This is handy for filling up volumetric flasks.

**Note 4:** At the beginning of the addition add 1 – 2mL  $\text{KMnO}_4$  solution. Then you have to wait a little and the deep violet color will disappear. After that you can continue with adding.

### Protocols

Write a protocol of the spectroscopic analysis on a blank sheet.

Write a protocol of your volumetric analysis on a blank sheet.


All found values must be noted clearly and in a reasonable order. Note your calculations in a reproducible way.

# 4-Methoxyacetophenone

## 8. Determination of 4-Methoxyacetophenone with External Standard on HPLC

A Friedel-Crafts reaction of anisole with acetyl chloride gives 4-methoxyacetophenone. In this case, an excess of anisole was used as solvent. The sample here is the raw product of 4-methoxyacetophenone with mainly anisole and a little tertiary butyl methyl ether from extraction as byproducts. It is your task to determine the amount of 4-methoxyacetophenone in the given sample.

## 9. Reaction Scheme of the Friedel-Crafts reaction



## 10. Reagents for HPLC Analysis

Reagent	Amount	Information	Safety advice
4-Methoxyacetophenone	ca. 2*50mg	Reference	H: 302-315 P: no advice
Raw 4-Methoxyacetophenone	ca. 2*100mg	Sample	<i>For Anisole:</i> H: 226 P: 210-262
Methanol HPLC grade	ca. 500mL	For dilution	H: 225-331-311-301-370 P: 210-233-280-302+352
Water nanopure HPLC grade	ca. 150mL	For dilution	No advice

## 11. Equipment for HPLC Analysis

Analytical balance	1 4all	Weighing tongs	1 4all
Volumetric flask with stopper 50mL	4	Piston stroke pipet 100μL	1 4all
Erlenmeyer flask with stopper 1000mL	1	Crimp capped sample vials	4
Measurement cylinder 250mL	1	Crimp cap tongs	1 4all
Spatula	2	Rack for vials	1
Piston stroke pipet 1000μL	1 4all		

Pasteur pipettes as required.

## 12. Procedure

### Preparation of diluent

In a stoppered 500mL Erlenmeyer flask, 200mL methanol *for HPLC* and 200mL water *nanopure* (for HPLC) are mixed until a clear solution is present.

### Dilution of the reference

Approximately 55mg 4-methoxyacetophenone are filled in a 50.0mL volumetric flask, weighed exactly and diluted with methanol / water nanopure 1:1 diluent to the mark. The solution is homogenized well. 625 $\mu$ L of the solution is pipetted with a 1000 $\mu$ L piston stroke pipette into a 25.0mL volumetric flask and diluted with methanol / water nanopure 1:1 diluent to the mark. This procedure is done twice and the flasks are reasonably labelled.

### Dilution of the sample

With the aid of a piston stroke pipet 100 $\mu$ L sample are dropped into a 50.0mL volumetric flask and the weigh is taken exactly. The sample is diluted with methanol / water nanopure 1:1 diluent to the mark. The solution is homogenized well. 625 $\mu$ L of the solution is pipetted with a 1000 $\mu$ L piston stroke pipette into a 25.0mL volumetric flask and diluted with methanol / water nanopure 1:1 diluent to the mark.

This procedure is done twice and the flasks are reasonably labelled.

### Preparation of vials

From each diluted volumetric flask a sample is taken and filled into a crimp vial. The crimp vial is tightly capped and clearly marked with a marker pen with sample1, sample2, reference1 and reference2. Then it is handed over to the HPLC operator. The chromatograms are taken by the HPLC operator with a prepared and tested method. You have to add them to the protocol as an appendix.

### HPLC parameters

Stationary Phase: NUCLEOSIL 120-5 C-18  
 Temperature: 35°C  
 Mobile phase: methanol / water 65% / 35% with 0.7mL/min flow, 1.2min stop time.  
 Injection: 10 $\mu$ L  
 Detection: 278nm

4-Methoxyacetophenon will give a clear, symmetric peak at about 0.7min with a height of 300 to 500mAU. Anisole will give a small symmetric peak at about 0.9min with a small height.

### Schedule for HPLC

One analysis run needs 2 minutes time from injection to injection. You have to hand over your samples before 2:30 p.m. for the HPLC analysis. Otherwise there might not be enough time for your samples to run.

### Calculations

Calculate the percentage of 4-methoxyacetophenon in the sample with the method of external standard; Once with reference1 and sample1 and once with reference2 and sample2.

Calculate the mean value of the two percentages.

Calculate the recovery with reference1 as reference and reference2 as a virtual sample. The found recovery percentage should be between 98 to 102%.

### Protocols and calculations HPLC analysis

Your protocol must contain a suitable table with the amounts of your sample and those of your references. In this table, write down the found areas and the results of your calculations. Describe your dilution work in a few words.

The calculation work can be done on a separate blank sheet or on a excel sheet on a computer.

Note your calculations in a reproducible way.

## 13. General Notes

You must pour all solutions containing silver in a special bin labeled "Ag<sup>+</sup> – waste". All other aqueous waste can be poured into the sink under dilution with water.

After use, rinse all glassware with demineralized (WBI) water. If there are chemicals visible, rinse it with ethanol first and with demineralized (WBI) water. If this procedure is insufficient, clean it with soap water and a test tube brush, and then rinse it well with tap water followed by demineralized (WBI) water.

All single use plastic devices are collected in the chemical waste box.